


SMART TAPER MICROLAPAROSCOPIC INSTRUMENT



SMART  TAPER

KEY FEATURES and BENEFITS

3-piece design for ease of cleaning



Threaded connections for easy and secure assembly

Smart Taper



Patented SMART TAPER instruments provide the jaw size and strength of 5mm instruments, through microlaparoscopic trocars.



SMART TAPER

SMART  TAPER

JAW INSERTS

Fenestrated Forceps

	33 cm length
20 mm	C35-IP-F20-33
30 mm	C35-IP-F30-33
30 mm CURVED	C35-IP-F3B-33
15 mm DOUBLE ACTION	C35-IP-2FS-33
20 mm DOUBLE ACTION	C35-IP-2F2-33
25 mm DOUBLE ACTION	C35-IP-2F3-33




Grasping Forceps

	33 cm length
CLUTCH	C35-IP-CL-33
CLUTCH 20 MM DOUBLE ACTION	C35-IP-2CL-33
CLUTCH 30 MM DOUBLE ACTION	C35-IP-2C3-33
STRAIGHT 15 MM	C35-IP-S15-33
BABCOCK	C35-IP-BK-33
INTESTINAL CLAMP	C35-IP-IC-33



Toothed Grasping Forceps

	33 cm length
ALLIS	C35-IP-AL-33
ALLIS DOUBLE ACTION	C35-IP-2AL-33
SEMM DOUBLE ACTION	C35-IP-2SM-33
BIOPSY DOUBLE ACTION	C35-IP-2BI-33



Dissectors

	33 cm length
MARYLAND DISSECTOR	C35-IP-2DM-33
S SHAPE DISSECTOR	C35-IP-2DS-33



Scissors

	Disposable sterile
CURVED SCISSORS TIP	C35-ST-2SC-ST



INSULATED OUTER TUBES and HANDLES

Insulated Outer Tubes

	33 cm length
FORCEPS INSULATED OUTER TUBE	C35-TU-33
SCISSORS INSULATED OUTER TUBE	C35-TUSC-33



Handles

	PEEK HANDLE WITHOUT RATCHET	PEEK HANDLE WITH RATCHET
YELLOW - GRASPING FORCEPS	LIL-H-YP1	LIL-H-YP2
BLUE - DISSECTORS	LIL-H-BP1	LIL-H-BP2
GREEN - TOOTHED GRASPING FORCEPS	LIL-H-GP1	LIL-H-GP2
RED - SCISSORS	LIL-H-RP1	/





» Color coding is to aid in identification; all handles are interchangeably

TROCARS and ACCESSORIES

Reusable Trocars with stopcock

	CANNULA	OBTURATOR
L=60mm	COEC35060RO	COEP35060H
L=75mm	COEC35075RO	COEP35075H
L=100mm	COEC35100RO	COEP35100H
L=120mm	COEC35120RO	COEP35120H



Disposable Sterile Valves

DISPOSABLE STERILE VALVES

Patented 3 point valve eliminates the need for secondary valve within cannula

BOX OF 25

C35-SEAL-1

Individually packed valves

25x



BOX OF 30

C35-SEAL-3

3 packed valves

10x



Scalpel Accessory Tool

Assures precise skin incision to minimize scarring as well as trocar slippage

C35-SAT-2



Use with 11 blade, to create optimal skin incision

NEEDLE HOLDER and SPECIALTY

Needle holder

33 cm length

C35-CT-33



Hook

33 cm length


C35-CH-33



Aspiration Cannula

33 cm length


C35-PC-33



Aspiration Needle

33 cm length

C35-PN-33



Coagulation Tip

33 cm length

C35-CT-33



LAPAROSCOPE and BIPOLAR


Microlaparoscopes

0 DEGREES

30 DEGREES

LIL-OPT-0

LIL-OPT-30



Bipolar forceps

INSERT

33 cm length


C35-IB-FS-33



OUTER TUBE

33 cm length

C35-BTU-33



HANDLE

33 cm length

C35-H-BP



CABLES

Monopolar cable

80100016

Bipolar cable

80100021



BIPOLAR CABLE

MONOPOLAR CABLE



ab medica s.a.s.

Les Petites Quarterées - 18100 Mery Sur Cher, France
tel. +33(0)248515993 - www.abmedica.fr

ab medica s.p.a.

Via J. F. Kennedy, 10/12 - 20023 Cerro Maggiore (Mi), Italy
tel. +39 02 93305.1 - www.abmedica.it